

ECE 105: Introduction to Electrical Engineering

Lecture 13

Bio 3

Yasser Khan

Rehan Kapadia

Pulse oximetry




ECG

Pulse Ox

Respiration

Blood Pressure


Temperature






$$SO_2 = \frac{C_{HbO_2}}{C_{HbO_2} + C_{Hb}}$$

- SO_2 The saturation of oxygen in blood,
- C_{HbO_2} Concentration of oxygenated hemoglobin (HbO_2),
- C_{Hb} Concentration of deoxygenated hemoglobin (Hb).




- Pulse oximetry measures blood oxygenation. Using spectrophotometry of absorptivity of blood at two distinct wavelengths, blood oxygen saturation is quantified.
- Can detect hypoxemia, ie. lower than normal blood oxygenation.





Oxy and deoxy-hemoglobin




Oxy and deoxy-hemoglobin




Transmission-mode oximetry






How the transmitted light is changing




Oxy and de-oxy abs.





Why do we need two lights




Change in ratio = change in oxygenation





Need for empirical calibration




Data collection flow




AC and DC part of the signal




AC and DC part




The graphs you will see in your demoboard




Pulsatile signal




The portion of the AC part is very small





Source of drift




Optical shunting – photodiode gets saturated








Transmission vs. reflectance oximetry




- Transmission-mode pulse oximetry is limited only to tissues that can be transilluminated, such as the earlobes and the fingers.
- If reflected light is used as the signal, the sensor can be used beyond the conventional sensing locations.







- AC signal is the highest at the forehead for both Red and NIR channels.
- Arms provide mid-range AC amplitude, while signal strength is low in the legs and chest area.
- Forehead is the best location for reflectance pulse oximetry.



Oximeter readout circuit




Calculating oxygen saturation





$$R = \frac{AC_{rd}/DC_{rd}}{AC_{ir}/DC_{ir}}$$

$$R = \frac{20mV/625mV}{20mV/330mV} = .528$$


Press capture




Keep your finger gently on the sensor




Going back to the PPG data collected




Smoothing




Drift correction




Peak detection




Heart rate




Using single component to find heart rate




4 components in the PPG signal – signal loss



Taking all the components provides the original signal



Ratio of the ratios



(Red PPG Signal) DC Average: 99986.91, AC Average: 139.15
(IR PPG Signal) DC Average: 116077.38, AC Average: 386.22
Red signal ratio (red_ac/red_dc): 0.0013916875411345753
IR signal ratio (ir_ac/ir_dc): 0.0033272520894769473
Ratio of Ratios (Red ratio/ IR ratio): 0.4182693416997304

